

The Analysis of 85 Pesticides Commonly Found in Cannabis and Cannabis-Related Products by LC-MS and GC-MS

By Abdullah Alsamraie, Kayla Dermody, Philip Brown, Susan Meronek

Abstract

Residual pesticides in cannabis and cannabis products can pose a major health risk to consumers in both recreational and medicinal applications. The main source of these chemicals comes during the growing process where chemicals are applied to the plants to control insects, mold or fungus growth. Some of those pesticides, such as organophosphates, are known for their toxic effects to both humans and wildlife.

As a result, regulators and standardization bodies such as ASTM and AOAC have published different test methods to analyze for residual pesticides in cannabis and cannabis-related products.

In order to make this screening as efficient as possible, AccuStandard used the target lists from different agencies, and was able to identify and quantify over 80 pesticides commonly used during cannabis cultivation. The analysis was conducted using multiple Certified Reference Materials (CRMs) prepared by AccuStandard and formulated based on the compatibility and the stability of those chemicals. Using both GC-MS and LC-MS we were able to offer a comprehensive and robust analysis for those compounds, most of which are in the target list of compounds in ASTM test method D8399 and AOAC SMPR 2018.011. The results illustrate the proper methods for handling the CRM samples and the appropriate analytical technique used to evaluate every compound. The conclusion is that using both analytical techniques combined with well-characterized CRMs represents a significant benefit for an accurate and reliable analysis of these harmful chemicals.

Introduction and Background

Pesticides are a group of chemical compounds that are primarily used to prevent and control harmful organisms and diseases during the plant growing process. Most of those chemicals are regulated by different agencies due to their harmful health effects. Cannabis is no different than other crops where pesticides are essential to improve quality and increase yield. While the USEPA and other federal agencies regulate pesticide residual limits on crops, no federal guidance has been established for pesticides and its residual content in cannabis due to its illegal federal status. Since residual pesticide analysis is a major part of ensuring cannabis is safe for human consumption, this has left consumers, especially individuals who use cannabis for therapeutic purposes, at risk of being exposed to a wide array of toxic chemicals. As a result, a reliable and accurate characterization of these contaminants is important for the safety of consumers and to reduce exposure.

In this study, AccuStandard was able to analyze and quantify 85 pesticides commonly found in cannabis and its related products. The analysis was conducted using 5 different Certified Reference Materials (CRMs) and using both LC-MS and GC-MS to examine method suitability. Compounds were accurately identified and confirmed using retention time and its respective mass spectrum.

This work also presents a comprehensive approach for the analysis of residual pesticides in cannabis to support labs in meeting the local and state regulatory guidelines. This has been accomplished by using suitable analysis techniques and reliable CRMs to help overcome the challenges of identifying these contaminants.

85 Targeted Compounds

Abamectin	Chlorantraniliprole	Fenoxycarb	Methiocarb	Resmethrin
Acephate	Clofentezine	Fenprophymate	Methomyl	Spinetoram
Acetamiprid	Clothianidin	Fensulfthion	Methoprene	Spinosad
Acetaminophol	Coumaphos	Fenitrothion	Mevinphos	Spirodiclofen
Alar	Cyrantraniliprole	Fenvalerate	Myclobutanil	Spiromesifen
Aldicarb	Cyprodinil	Fipronil	Naled	Spirotetramat
Allethrin	Deltamethrin	Flonicamid	Novaluron	Spiroxamine
Azadirachtin	Diazinon	Fludioxonil	Oxamyl	Sumithrin
Azoxystrobin	Dichlorvos	Fluopyram	Pacllobutrazol	Tebuconazol
Baygon	Dimethoate	Hexythiazox	Permethrin	Tebufluozide
Benzovindiflupyr	Dimethomorph	Imazalil	Piperonyl butoxide	Teflubenzuron
Bifenazate	Dinotefuran	Imidacloprid	Pirimicarb	Tetrachlorvinphos Z
Bifenthrin	Dodemorph acetate	Imidan	Prallethrin	Tetramethrin
Boscalid	Dursban	Iprodione	Propiconazole	Thiacloprid
Buprofezin	Ethopropr	Kresoxim-methyl	Pyraclostrobin	Thiamethoxam
Carbaryl	Etofenprox	Malathion	Pyrethrins	Thiophanate-methyl
Carbofuran	Etoazole	Metaxyl	Pyridaben	Trifloxystrobin




CRM Handling and Methods of Analysis

CRMs were designed and prepared to maximize stability, compatibility, and shelf life. Some compounds such as Spirodiclofen and Spiromesifen were prepared independently from the multiple component mixes due to their unstable nature. Mixes were sonicated for 45 minutes and acclimated to room temperature prior to analyzing to ensure solubility and improve each component's response. The concentration of all components was verified against another independently made standard.



LC Parameters

- General gradient method starting @ 100% water for 1 min @ 0.3 mL/min, the 100% Acetonitrile in 13 min, 100% Acetonitrile until T=16 min
- Column Zorbax Eclipse plus C-18, 2.1 mm X 50 mm, 1.8 µm particle size
- Agilent 1290 Infinity II LC/MS system
- Agilent 6125 quadruple MS detector




Peak # Component	Peak # Component	Peak # Component	Peak # Component	Peak # Component
6.075 Imidacloprid (138261-41-3)	4.849 Oxamyl (23135-22-0)	0.593 Alar (1596-84-5)*	4.496 Dinotefuran (165252-70-0)	13.389 Spirodiclofen (148477-71-8)
7.159 Imazalil (35554-44-0)	4.980 Methomyl (16752-77-5)	3.621 Acephate (30560-19-1)*	5.394 Pirimicarb (23103-98-2)	
6.690 Spinosad (168316-95-8)	5.392 Flonicamid (158062-67-0)	6.188 Dimethoate (60-51-5)	5.883 Clothianidin (210880-92-5)	
9.029 Spinetoram (168316-95-8)	5.475 Thiamethoxam (153719-23-4)	6.936 Aldicarb (116-06-3)	6.448 Mevinphos (7786-34-7)	
9.326 Spirotetramat (203313-25-1)	6.323 Acetamiprid (135410-20-7)	6.936 Thiacloprid (111988-49-9)	7.403 Dichlorvos (62-73-7)	
9.604 Myclobutanil (88671-89-0)	8.158 Carbaryl (63-25-2)	7.588 Metatol (31717-87-0)*	7.723 Thiophanate-methyl (23564-05-8)	
9.733 Azoxystrobin (131860-33-8)	8.278 Metalaxyl (57837-19-1)	7.743 Baygon (114-26-1)	8.451 Cyrantraniliprole (736994-63-1)	
9.820 Tebuconazol (107534-96-3)	8.588 Clofentezine (74115-24-5)	7.856 Carbofuran (1563-66-2)	8.451 Cyprodinil (12152-61-2)	
10.078 Bifenazate (149877-41-8)	8.906 Dimethomorph (110488-70-5)	7.856 Spiroxamine (118134-30-8)	8.544 Fensulfthion (115-90-2)	
10.278 Malathion (121-75-5)	9.077 Chlorantraniliprole (500008-45-7)	8.129 Azadirachtin (11141-17-6)	9.385 Fluopyram (658066-35-4)	
12.665 Etoazole (153233-91-1)	9.027 Dimethomorph (110488-70-5)	8.785 Naled (300-76-5)	10.069 Buprofezin (69327-76-0)	
13.092 Abamectin (17151-41-2)	9.147 Spinetoram (935545-74-7)	9.047 Pacllobutrazol (76738-62-0)	10.202 Iprodione (36734-19-7)	
13.692 Permethrin (52645-53-1)	9.471 Spinetoram (935545-74-7)	9.276 Methiocarb (2032-65-7)	10.223 Tetrachlorvinphos Z (22248-79-9)	
13.930 Permethrin (52645-53-1)	9.512 Fludioxonil (131341-86-1)	9.686 Ethopropr (13194-46-4)	10.563 Tebufenozide (112410-23-8)	
	9.804 Boscalid (188425-85-6)	9.686 Imidan (732-11-6)	11.005 Fenitrothion (55-38-9)	
	10.279 Tilt (60207-90-1)	10.282 Fenoxycarb (72490-01-8)	11.213 Pyraclostrobin (175013-18-0)	
	10.775 Kresoxim-methyl (143390-89-0)	10.794 Fipronil (120068-37-3)	11.267 Coumaphos (56-72-4)	
	11.233 Clofentezine (74115-24-5)	10.794 Benzovindiflupyr (1072957-71-1)	11.444 Allethrin (584-79-2)	
	11.682 Trifloxystrobin (141517-21-7)	11.046 Diazinon (333-41-5)	11.632 Novaluron (115714-46-6)	
	12.492 Hexythiazox (78587-05-0)	11.531 Pyrethrins (8003-34-7)	12.050 Tetramethrin (7696-12-0)	
	12.622 Fenprophymate (111812-58-9)	12.160 Pyrethrins (8003-34-7)	12.160 Tetramethrin (7696-12-0)	
	12.046 Piperonyl butoxide (51-03-6)	12.464 Dursban (2921-88-2)	13.453 Resmethrin (18459-86-8)	
	12.464 Dursban (2921-88-2)	12.917 Pyrethrins (8003-34-7)	13.538 Deltamethrin (52918-63-5)	
	13.439 Pyrethrins (8003-34-7)	13.439 Pyrethrins (8003-34-7)	13.538 Fenvalerate (51630-58-1)	
	13.950 Etofenprox (80844-07-1)	13.950 Etofenprox (80844-07-1)	14.050 Methoprene (40596-69-8)	
	14.260 Bifenthrin (82657-04-3)	13.950 Etofenprox (80844-07-1)		
	15.399 Acequinolyl (157960-19-7)	15.399 Acequinolyl (157960-19-7)		

* Compounds detected on MS

GC Parameters

- Cons. Split, 12psi constant flow: Split 20:1, 1 µL inj; GCMS, INJ 270C
- Column: DB-5MS 30m, 0.25iD, 0.25 µm
- Oven program: 80°C, 17°C/min to 340°C, 8min



Peak # Component	Peak # Component	Peak # Component	Peak # Component	Peak # Component
2.210	3.344	4.947	6.538	7.984
8.042	8.255	9.367	9.710	10.426
10.267	10.426	10.426	10.426	10.426
11.281	11.281	11.281	11.281	11.281
11.904	11.904	11.904	11.904	11.904
12.571	12.571	12.571	12.571	12.571
13.571	13.571	13.571	13.571	13.571
14.271	14.271	14.271	14.271	14.271
14.433	14.433	14.433	14.433	14.433
15.022	15.022	15.022	15.022	15.022

Additional Notes

- Standards were designed based on the stability and compatibility of each individual compound
- Each mix was prepared in Acetonitrile at 100 µg/mL
- All analytes were confirmed using Mass Spectrometry



Conclusion

As more states legalize cannabis for both medicinal and recreational uses, cannabis testing labs face growing pressure to test samples for more pesticides efficiently and accurately to ensure the safety of consumers. During this work we were able to demonstrate the effectiveness of using appropriate and accessible analytical methodologies, in addition to well-characterized Certified Reference Materials (CRMs), in the identification and quantification of residual pesticides in cannabis. Using a combination of both LC-MS and GC-MS techniques offered a significant advantage for the selection of suitable methods to reliably and accurately analyze all targeted compounds in this study.